

## Track-N

Features: Non-Compliance Balloon Catheter

- \* Tailor- made polymer with sleek finish to refrain collision with stent struts
- \* Direct torque transmission design allows superior steerability and control
- \* Pushability, trackability and excellent progressive flexibility of the distal area
- \* Ultra smooth transition area and lubricious hydrophilic coating for friction free performance
- \* Ergonomic, transparent hub facilitates safe connection
- \* Proximal marker bands at 90 and 100cm
- \* A parallel guide wire exit port guarantees best pushability and a better dye flow injection
- \* Lower profile, new jacketed hypo tube design and an exclusive blend of distal shaft material for more responsive push and track
- \* Excellent, responsive trackability by virtue of its design
- \* Highly flexible catheter outer body maximizes trackability performance
- \* Focused inflation force
- \* Controlled pushing of stent towards the vessel wall
- \* Incomplete apposition may contribute to thrombosis formation and SAT
- \* Post Dilation reduces target vessel revascularization (TVR)
- \* Uniform stent apposition facilitates uniform drug absorption into endothelial tissue

## Track-N Non-Compliance Balloon Catheter

Technical Specifications							
Design	Rapid Exchange						
Balloon compliance	Non- compliant						
Distal Shaft (outer diameter)	0.90 mm						
Proximal Shaft (outer diameter)	0.72 mm						
Usable catheter length	1414 mm						
Crossing profile	0.8 - 1.4 mm						
Minimum guide catheter compatible	0.058" / 5 Fr						
Guide wire compatible	0.014'						
Nominal inflation pressure	10 ATM						
Rated Burst Pressure	18 ATM						
Marker bands	Two platinum iridium markers						
Tip profile	0.016"						

Inflation Pressure (ATM)	Balloon Diameter by Size (mm)														
	1.50	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00	4.25	4.50	4.75	5.00
6	1.46	1.70	1.94	2.19	2.41	2.68	2.93	3.18	3.41	3.68	3.90	4.17	4.40	4.66	4.90
8	1.48	1.73	1.97	2.21	2.45	2.72	2.96	3.21	3.45	3.72	3.95	4.21	4.45	4.70	495
10*	1.50	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00	4.25	4.50	4.75	5.00
12	1.52	1.77	2.03	2.28	2.54	2.79	3.05	3.30	3.55	3.79	4.05	4.30	4.56	4.80	5.08
14	1.55	1.79	2.06	2.32	2.57	2.83	3.10	3.34	3.60	3.84	4.11	4.35	4.61	4.86	5.15
16	1.57	1.82	2.09	2.35	2.61	2.87	3.15	3.37	3.65	3.89	4.16	4.40	4.67	4.92	5.22
18**	1.59	1.85	2.12	2.39	2.65	2.92	3.20	3.41	3.69	3.94	4.20	4.45	4.73	4.98	5.29
20	1.62	1.88	2.16	2.42	2.68	2.97	3.25	3.44	3.73	3.99	4.25	4.50	4.79	5.03	5.36
22	1.65	1.91	2.20	2.45	2.72	3.01	3.29	3.48	3.78	4.03	4.29	4.54	4.85	5.09	5.43

■ Nominal Pressure ■ Rated Burst Pressure (RBP)



## Sah Health Care Pvt. Ltd.

POCKET C- 565 PAPER MARKET IFC, GHAZIPUR MAYUR VIHAR PHASE - 3 DELHI 110096 Mob No +91 9990045799 Email id sahhealthcarepvtltd@gmail.com