



SAH TRACK™

Track-T Wire

Hydrophilic Guidewire



Product Characteristic

- Acqua wire has non traumatic, flexible tapered tip which increase flexibility for excellent smooth and safe navigation through vessels.
- Improved hydrophilic coating provide smoother navigation through catheter and vessels. It also prevents blood platelet adhesion and thrombus formation.
- With new polyurethane radiopaque jacket Acqua wire has very smooth surface to minimize blood adhesion to the guidewire.
- Rounded end decrease likelihood of vessel traumas smoother guide wire insertion.

Ordering Information

Standard Hydrophilic Guidewire : 150cm

REF	Specification	Coating
NT/GWH/32150ST	0.032 Inch, 150cm, ST Tip	Hydrophilic Coating
NT/GWH/35150ST	0.035 Inch, 150cm, ST Tip	Hydrophilic Coating
NT/GWH/32150A	0.032 Inch, 150cm, Angle Tip	Hydrophilic Coating
NT/GWH/35150A	0.035 Inch, 150cm, Angle Tip	Hydrophilic Coating

*Standard Configuration : 3mm St/Angle Tip, Standard Stiffness, Fixed Core, Hydrophilic Coating

Standard Hydrophilic Guidewire : Small Size

REF	Specification	Coating
NT/GWH/32260ST	0.032 Inch, 260cm, ST Tip	Hydrophilic Coating
NT/GWH/35260ST	0.035 Inch, 260cm, ST Tip	Hydrophilic Coating
NT/GWH/32260A	0.032 Inch, 260cm, Angle Tip	Hydrophilic Coating
NT/GWH/35260A	0.035 Inch, 260cm, Angle Tip	Hydrophilic Coating

*Standard Configuration : 3mm St/Angle Tip, Standard Stiffness, Fixed Core, Hydrophilic Coating

Hydrophilic Guidewire

Track-T Wire is a Standard type are Nitinol hydrophilic guidewires covered with polyurethane and a hydrophilic polymer coating for an excellent manoeuvrability and absence of friction. Indicated for standard angiography in Coronary and Peripheral territories, catheterization, embolization, peripheral angioplasty and stenting (with or without exchange wire), and subintimal angioplasty in small, standard or large proximal or distal vessels.

Sah Health Care Pvt. Ltd.

Mob No +91 9990045799 Email id sahhealthcarepvtltd@gmail.com